Instruction for use Physics® Forceps

Intended use:

Tooth extraction. Separate the gingival attachment from the tooth.

Contraindication:

Impacted wisdom tooth, ankylosed tooth and tooth with non-solid root.

Prior to initial use:

Always read thoroughly the instructions for use.

Improper use such as misplaced bumper and excessive force can cause root-tip fracture.

Webpage: www.directadental.com/luxator/physicsforceps

Precaution


- Inspect product for proper function before use. Broken product can cause injuries and shall not be used.
- Inspect surface of product for corrosion. Corroded product can cause cross-contamination and shall not be used.
- Re-usable instruments shall be processed according to processing instruction in between uses to prohibit cross-contamination
- Always use bumper guards to avoid soft tissue injuries
- Lower molar (tooth with multiple root tips) is recommended to be sectioned before removal, to ease removal and avoid root tip fracture

Warning


- Do not re-use bumpers, those are for single-use as they cannot be cleaned properly and may have damages after use.
- Misplaced beak gives less lever effect
- Misplaced bumper can cause less lever effect and more force will be required, which can cause root-tip fracture, the force shall be perpendicular to the root
- Do not squeeze the handles as this can cause damage of alveolar bone

How to use


Step 1


Step 2

Step 3

Step 4

Step 5

physics®
forceps

EN

Physics® Forceps Extraction instruments

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DIRECTA INC
64 Barnabas Road, Unit 3
Newtown, CT 06470, United States
203.491.2273

Rx only

CE Made in USA

DIRECTA

106715-2022-09-20

Please note that this step is not required.

Separate the gingival attachment from the tooth. No elevation is required prior to using the Physics Forceps. Consider use of Luxator Periotope or Luxator P-series in preparation for extraction.

Step 1

This step is not always necessary, but it helps when the tooth is decayed for a better purchase point. Use a flame shaped diamond bur to drill a trench on the lingual side of the tooth (at the expense of the tooth) just below the gum line to secure a solid purchase point with the beak. If utilizing the Molar Series and placing the bumper on the lingual side, this trench would then be on the buccal side of the tooth.

Step 2

With the handles opened wide, set the beak into the depth of the lingual sulcus (or optional buccal sulcus with Molar Series) as far as possible on solid root surface. In general, the beak is placed before the bumper.

Step 3

Set the bumper perpendicular to the tooth at about the level of the mucogingival junction, or as deep in the vestibule as possible. In utilizing the Molar Series with the bumper on the lingual, place the bumper on the palate for the uppers or as deep in the vestibule as possible for lowers. Hold that position securely. **STOP! DO NOT SQUEEZE** the handles.

Step 4

Without squeezing the handles or moving your arm, begin to apply a steady, very slow rotational force towards the buccal (or optionally toward the lingual with Molar Series). As soon as you get resistance, STOP and hold the instrument in that position for **15-20 seconds**. After this time if the tooth is not moving, move your wrist a little further and hold again for **15-20 seconds**. Continue this pattern. Be patient as this can take up to 2-3 minutes with no movement. The force, or creep, will continue to build, allowing the hyaluronic acid to build and causing the PDL to release.

Step 5

As soon as the tooth lifts up and out of the socket 1-3 mm, the Physics Forceps has completed its designed function. Deliver the tooth with Luxator Forte or a hemostat, rongeurs or conventional forceps. Do NOT rotate the tooth all of the way out of the socket when using the Physics Forceps.

Processing instructions for re-usable instrument only

Warnings	<p>Improper cleaning/disinfection and/or sterilization can corrode, pit, crack or stain instruments. Generally, a pH less than six is acidic and may breakdown the protective surfaces of stainless steels, resulting in pitting and/or black staining; and a pH more than eight is alkaline and can cause brown staining which can also interfere with the smooth operation of instruments. (Most brown stains are not rust and are easily removed with surgical stain remover.)</p> <p>Do not use steel wool, drill brushes or abrasive items.</p> <p>Do not batch stainless steel, aluminium, brass or copper instrument together. If batched together, a potential for electrolysis reactions between dissimilar metals will exist, which can produce etching and corrosion on the instrument surface.</p> <p>Instrument should not be in contact with the following chemicals for more than a couple of hours (then immediately and thoroughly rinsed): aluminium chloride, barium chloride, mercury dichloride, calcium chloride, carbolic acid, citric acid, cresol, mercury chloride, mercury salts, phenol, permanganic acid potash, potassium thicyanate, ferrous chloride, stanniferous chloride, tartaric acid.</p> <p>The following chemicals should be avoided completely: Aqua Regia, iron chloride, sulphuric and hydrochloric acid and iodine.</p>
Limitations on reprocessing	See inspections.

Instructions

Initial treatment at the point of use	Remove the Bumper Guard from the instrument.	
Preparation before cleaning	Instruction	Wipe off the instrument after use to prevent soil and debris to dry onto the instrument. Perform cleaning as soon as possible after use. Do not exceed 2 hours.
Cleaning: Manual	Equipment	Soft bristle brushes of various sizes Ultrasonic bath
	Detergent	Enzymatic or low alkaline, suited for medical devices. Follow agent manufacturers recommendations for concentration and temperature. Ultrasonic bath: non-corrosive, neutral cleaning agent with minimal foaming characteristics.
	Water quality	Drinking water quality for manual cleaning and distilled or de-mineralized water for ultrasonic bath and final rinse should be used.
	Instruction	<ol style="list-style-type: none"> 1. Immerse in freshly prepared cleaning solution as per manufacturer's instructions although maximum 40°C. Clean mechanically with a brush, working beneath the liquid level, until visibly clean. 2. Thorough rinsing 3. Clean in an ultrasonic bath for a minimum of 5 minutes, using a frequency of 35-55 KHz and power min. 50 W. Make sure the instruments are in an open position. 4. Rinse thoroughly for min 30 seconds. 5. Carefully dry instruments with clean compressed air (Class 1 or better according ISO 8573-1:2010) or lint free wipe. Any remaining water can cause potential rust or corrosion on the instrument surface and into the substrate. 6. Inspect cleanliness. If necessary, repeat from step 1.
Cleaning: Automated	Equipment	Washer-disinfector
	Detergent	Enzymatic or low alkaline, suited for medical devices. Follow detergent manufacturer's instructions regarding concentration and temperature.
	Rinsing Agent	Suited for medical devices. Follow agent manufacturer's instructions regarding concentration and temperature.
	Water quality	Drinking water quality for cleaning and deionized or purified water for final rinse/disinfection.
	Instructions	<ol style="list-style-type: none"> 1. Pre-wash the instruments under cold running water, use a soft bristle brush to remove visible dirt, pay extra attention to the joints and areas hard to access. 2. Load the instruments in an open position in the washer-disinfector. 3. Run program: <ul style="list-style-type: none"> - Pre-wash in cold water, 5min - Main wash 65°C, 10 min - Rinse in warm water 55°C -65°C, 3 min. - Rinse in warm water and rinsing agent 55°C-65°C, 10 min - Final rinse in DI or PURW water 65°C, 5min. - Drying at 110° C, minimum 15 min 4. When unloading, control that the instruments are clean. If necessary, repeat from step 1 or use manual cleaning
Disinfection: Manual	Detergent	Immerse into disinfection solution suited for instruments of stainless steel. Follow the disinfection solution manufacturer's instruction for concentration and time. If necessary, rinse instruments thoroughly with RO or distilled water, to remove remaining disinfection from the surface. Place the instruments on a tray for drying.

Disinfection: Automated	Equipment Water quality Detergent	Washer-disinfector (EN ISO 15883) Drinking water quality for cleaning and deionized or purified water for final rinse/disinfection. Enzymatic or low alkaline, suited for medical devices. Non-corrosive, neutral rinsing agent with minimal foaming characteristics. Follow agent manufacturer's instructions regarding concentration and temperature.
	Instruction	<ol style="list-style-type: none"> 1. Pre-wash the instruments under cold running water, use a soft bristle brush to remove visible dirt, pay extra attention to the joints and areas hard to access. 2. Load the instruments in an open position in the washer-disinfector. 3. Run program: <ul style="list-style-type: none"> - Pre-wash in cold water, 5 min - Main wash 65°C, 10 min - Rinse in warm water 55°C -65°C, 3 min - Rinse in warm water and rinsing agent 55°C -65°C, 10 min - Final rinse in DI or PURW 93°C, 5 min - Drying at 110° C, minimum 15 min 4. When unloading, control that the instruments are clean. If necessary, repeat the cycle or use manual cleaning.
Inspection, maintenance and lubrication	Inspect surface of product for corrosion or other damages. Corroded product can cause cross-contamination and shall not be used. Lubricate the joint with approved lubricants intended for medical devices.	
Packaging for sterilization	Equipment Instruction	Standard packaging pouch Pack the instrument in an open position in a sterile bag
Sterilization	Caution Equipment	The instrument must be cleaned and disinfected before sterilization. Steam autoclave (EN 13060 or EN 285 or EN ISO 17665-1)
	Instructions	Run minimal cycle: Steam temperature Min. 134 °C / 273°F Steam exposure time 4 min Vacuum drying 3 min
Storage	Ensure the instruments are dry before storage. The instrument should be stored in a packed tray, free from dust and moisture at moderate temperature of 5°C – 40°C	
Transportation	Not applicable	
Additional information	All equipment, detergents, lubricants and disinfections solutions shall be CE marked as per 93/42/EEC or REGULATION (EU) 2017/745	

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DIRECTA AB