ENGLISH

Disposable dental splint for hospital use

Using blunt instruments in the oral cavity in connection with peroral endoscopy or intubation can be hazardous to the teeth. This occurs under excitation or when instruments are pressed against teeth. On some occasions exfoliated teeth have entered the trachea or even a lung, thus making open chest surgery necessary (Acta Otolaryngol 84:450-52. 1977, McCarthy and Carlson). One way to minimize such risks is to use a dental splint.

CAMO disposable dental splints are made of a semi-hard external shell filled with soft mouldable silicone. This reinforced silicone base prevents the teeth from penetrating the mould. The integrity of the silicone is ensured by a thin plastic film, which prevents silicone fragments from entering the oral cavity. The size of the splint has been constructed to fit 90% of all jaw sizes. The size can, however, be easily adjusted with scissors.

- CAMO disposable splints can be used both in the upper and lower jaws.
- Both jaws can be protected simultaneously without excessively diminishing the necessary area for peroral endoscopy.
- The splints absorb high pressure without damaging teeth, thus the splints may be used during all kinds of pathological or induced excitations.

Note

Single use only. Shall not be sterilized or re-used as intended use can be jeopardized.

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