

DIRECTA ceramir[®] cement

THE NATURAL CHOICE FOR LASTING CEMENTATION

Easy to Use. Reliable. Bioactive.



Bioactive Strength For Lasting Restorations

Ceramir Crown & Bridge is a self-setting bio-compatible, and bioactive dental cement for conventional permanent cementation in various types of restorations.

Its powder consists mainly of calcium aluminate, and the liquid is water. The cement hardens through a chemical reaction. Once set, it reacts with saliva to form hydroxyapatite at the margins, creating a tight, permanent seal that minimizes microleakage and marginal gaps.

Its alkaline nature and natural integration enhance tissue compatibility, reduce post operative sensitivity, and lower the risk of secondary decay, resulting in a long-lasting patient friendly restoration. Ceramir dental cement requires no additional preparation of the tooth or the restoration before cementation.

It has a smooth flowability and generous working time, ensuring easy placement. During setting, it enters a controlled gel phase, allowing easy removal of excess cement. These properties make Ceramir dental cement an excellent choice for durable natural restorations.

Ceramir Crown & Bridge is offered in two forms: the pre-mixed, ready-to-use QuikCap capsule, and the QuikMix option, which allows you to mix the powder and liquid manually.

Indications


Indicated for permanent cementation of prosthetic restorations including metal, PFM (porcelain fused to metal), Zirconia and other high strength ceramics.



Biocompatible &
Bioactive




Generous working
time



Bioactive permanent
seal



Easy removal of
excess cement




No additional
preparation

"Ceramir Crown & Bridge has been the obvious choice at our clinic. It is a very easy-to-use cement with simple removal of excess material, while allowing the patient to bite down during the initial setting phase. Unlike competitors, one capsule is more than enough for a three-unit bridge – which is both convenient and saves money. We really recommend all our colleagues to use Ceramir Crown & Bridge."

- Tom Sandelin, Dentist, Kristallen Dental Clinic, Uppsala Sweden

Easy Handling for Effortless Placement



Easy to use, reliable results


Its smooth, flowable consistency and generous working time allow for effortless placement of the restoration. No additional preparation of the restoration or the tooth is required, saving you time and money. Ceramir dental cement is hydrophilic, wets and flows well even in moist conditions to ensure a smooth placement and adaptation of the restoration. It provides a strong, long-lasting

seal at the restoration interface, minimising the risk of secondary decay and offering long-term protection. Excess cement can easily be removed thanks to the controlled gel phase. These excellent handling properties minimise technique sensitivity and save valuable chair time.



Excess cement can easily be removed thanks to controlled gel phase.

Superior Sealing for Long-lasting Result




A permanent seal for lasting results

Ceramir dental cement provides a permanent seal at the restoration interface, reducing the risk of secondary decay and ensuring long-term protection. By eliminating gaps between the tooth and the restoration, it prevents bacterial leakage and supports the restoration over time, thereby contributing to its



longevity and a positive prognosis. Thanks to its alkaline properties (pH > 7), the cement promotes remineralisation and hydroxyapatite formation, strengthening the sealed interface and supporting the tooth's natural repair processes.

Ceramir Crown & Bridge	Baseline	7 days	23 days
Self-adhesive resin cement	Baseline	7 days	31 days
Resin-modified glass ionomer cement	Baseline	6 days	31 days
Glass ionomer cement	Baseline	11 days	58 days



Ceramir Crown & Bridge
Dentin

Ceramir cement Bond strength



Ceramir cement bond strength to different substrates. Data on file.

Reliable and durable adhesion

Achieving reliable and durable adhesion is essential for long-term clinical success. Ceramir dental cement delivers consistent bond strength comparable to leading resin-modified glass ionomer cements, ensuring secure retention across a variety of restorative materials.

Biocompatible & Hydrophilic - Inspired by Nature



Long-lasting, natural restorations

Ceramir dental cement supports biomimetic dentistry because its formulation is designed to work in harmony with the tooth. It is non-acidic, which makes it gentle on teeth and eliminates the risk of postoperative sensitivity. Free from resin, the cement is safe and stable, ensuring there is no shrinkage. Its natural bioactivity makes it biocompatible and hydrophilic, promoting remineralisation and forming hydroxyapatite to create a strong seal and long-lasting protection.

Ceramir cement in comparison

Property	IDEAL Material	Glass Ionomer	Resin modified Glass ionomer	Phosphate cement	Ceramir Cement
Film thickness µm	Low < 40	13,5	23,1	25,3	12-13 or < 25
Compressive Strength Mpa	High	122-162	40-141	62-101	150-200
Pulp irritation	Low	High	High	Moderate	Very Low
Removal of excess	Easy	Medium	Medium	Easy	Extremely Easy

Data on file.

Ceramir Bioceramic Implant Cement

Easy handling for implant treatments

Ceramir Bioceramic Implant Cement was specially developed to help dentists focused on implant treatments. Ceramir Bioceramic Implant Cement is a permanent bioceramic cement, 100% resin-free and tissue-friendly with excellent handling properties. No additional preparation of the

implant or restoration is required, saving you time and money. Ceramir Bioceramic Implant Cement ensures a smooth placement and adaptation of the restoration. Excess cement can easily be removed thanks to the controlled gel phase.




The main component in the powder is calcium aluminate and the liquid, water. The cement, which hardens because of a chemical reaction, is self-setting. In Ceramir Bioceramic Implant Cement QuikCap, the cement is supplied in a capsule for mixing and direct application. A capsule applicator and a mixing device are required for use. We recommend Ceramir Applicator 2.

Ceramir Applicator 2

A premium applicator designed for use with Ceramir QuikCap capsules for Ceramir Crown & Bridge and Ceramir Bioceramic Implant Cement.



Easy handling in a few steps



- 1. Activation**
Hold the capsule firmly against a hard, flat surface with the plunger facing down. Push down firmly, forcing the plunger into the body of the capsule until it stops.
- 2. Mixing**
Immediately place the capsule in a high-frequency oscillating (4,000 to 5,000 rpm) dental triturator or laboratory centrifuge. Mix for 8 seconds.
- 3. / 4. Application**
Immediately insert the activated capsule into the applicator and apply the cement. After application, remove the empty capsule. Working time (from end of mixing) at 23 °C/73 °F: 2 minutes.
- 5. Seating**
Seat the restoration and stabilize it under pressure, to avoid any movements, until the cement reaches a rubber-like consistency (approx. 2 min).
- 6. Removal of excess cement**
Remove excess material when the cement has reached the rubber-like consistency.
- 7. Finishing**
Keep the restoration stabilized for an additional 4 minutes. Check the area around the margins of the restoration (especially the gingival sulcus) and remove any remaining cement.



Design By Dentists

Directa's products are developed, tested and evaluated by our selected team of highly-qualified and renowned Swedish and international dental professionals. The aim is to make life easier for all dental professionals in their daily clinical work.

DIRECTA



DIRECTA AB

Finvids väg 8
194 47 Upplands Väsby, SWEDEN
+46 8 506 505 75
info@directadental.com
www.directadental.com

DIRECTA INC

64 Barnabas Road, Unit 3
Newtown, CT. 06470, USA
+1 203.491.2273
infousa@directadental.com
www.directausa.com



directadental



directadental



DirectaAB

The content of this brochure, including layout and images, is protected by copyright. ©2025 Directa AB. All rights reserved.

REF	Product Description	Pieces/Content
40030	Ceramir® Crown & Bridge QuikCap	20 capsules á 0,17 ml
40029	Ceramir® Crown & Bridge QuikCap	5 capsules á 0,17 ml
40032	Ceramir® Crown & Bridge QuikMix Kit*	Complete Kit
40033	Ceramir® Crown & Bridge QuikMix Powder refill	1 bottle of powder 15 g
40034	Ceramir® Crown & Bridge QuikMix Liquid refill	1 bottle of liquid 9 ml
40037	Ceramir® Bioceramic Implant Cement QuikCap	10 capsules á 0,17 ml
40027	Ceramir® Applicator 2	1 pcs/box

*1 bottle of powder 15 g, 1 bottle of liquid 9 ml, 1 mixing pad, 1 spoon

Reference sources:

A Review of Luting Agents

Cornelis H. Pameijer, International Journal of Dentistry Volume 2012, Article ID 752861, 7 pages
doi:10.1155/2012/752861

A bioactive dental luting cement- its retentive proper and 3- year clinical findings

Steven R. Jefferies, MS, DDS, PhD; Cornelis H. Pameijer, DMD, DSc, PhD; David C. Appleby, DMD, MScD, FACP; Daniel Boston, DMD; and Jesper Lööf, PhD Conexiant dental journal February 1, 2013 Volume 34 , Issue 2 , February 2013

Long-term retention of lithium disilicate crowns with a current bioactive cement

Kerry R. Streiff DDS | Xavier Lepe DDS, MS | Glen H. Johnson DDS, MS Wiley Received: 4 December 2020
Accepted: 19 January 2021 DOI: 10.1111/jerd.1271

Hydroxyapatite formation on a novel dental cement in human saliva

Johanna Engstrand 1, Erik Unosson, Håkan Engqvist
ISRN Dent. 2012:2012:624056. doi: 10.5402/2012/624056. Epub 2012 Sep 27.

Preliminary Evidence That Bioactive Cements Occlude Artificial Marginal Gaps

STEVEN R. JEFFERIES, MS, DDS, PhD*, ALEXANDER E. FULLER, DMD†‡, DANIEL W. BOSTON, DMD§
Journal of Esthetic and Restorative Dentistry, 2015 Wiley Periodicals, Inc. DOI 10.1111/jerd.12133



Ceramir dental
cement