Sendoline S5 Retreatment Files

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(en)

Intended purpose

Sendoline instruments covered by these instructions are intended for the removal of obturation materials from the root canal by dentists/endodontists in a clinical environment.

Design

The device operative part is made from nickel-titanium and is defined as a type 2 instrument according to EN ISO 3630-1. The shank is type 1 according to EN ISO 1797.

Contra-indications

General contra-indications for root canal treatment. Previously root canal treated teeth showing a severe curvature of the anatomical shape of the root canal.

Precautions

– The instruments are intended for single use with the purpose to avoid cross-infection and to minimize the risk of file fracture following repeated use.

– The instruments are sharp and must be handled with care to avoid injuries. Follow national healthcare regulations for safe disposal of sharps.

In case of a serious incident from using the instruments, immediately contact Directa AB and the applicable authority in your country.

Indications for use

The instruments are intended for mechanical instrumentation and should be used at the torque values (Ncm) and speeds (rpm) specified in the table found in this document.

– Create straight-line access to the canal orifice(s).

– The instruments shall be used in accordance with the recommended sequence, starting with file no 1 finishing with file no 4.

Sterilization instructions

Single use

The instruments are intended to be processed (sterilized) once prior to use and must not be reprocessed (re-sterilized) or reused.

Packaging

Arrange the instruments in an instrument tray or stand. Use suitable sterilization wrapping material or sterilization pouches compliant with EN ISO 11607-1. The packaging material must be large enough to ensure that the seal is not under tension. Do not overload the sterilization equipment.

Sterilization process

– Steam sterilizer, class B autoclave.
Applicable standards: EN 13060, EN 285, EN ISO 17665.

– Sterilization cycle type B in accordance with EN 13060.

– Sterilization holding time: minimum 3 minutes at 134°C (273°F).

– Vacuum drying: minimum 2 minutes at 134°C (273°F).

Do not exceed 136°C (277°F) during sterilization. The operator assumes responsibility for the use of any other sterilization procedures

Storage

Keep the devices in their sterilization packaging in a dry, dust free and clean environment.

Symbols

R_xonly For professional use only (USA)

 File type S - the S figurative mark is the property of Directa

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© 2025 Directa AB. All rights reserved.