


FenderMate® D-Coat

Matrix for Class II composite fillings

DIRECTA
directadental.com


1




2




3



4



5



FenderMate® D-Coat

left regular

left narrow

right regular

right narrow

t= 0,05mm

Food approved plastic (Polypropene) and

colours. Stainless steel.

DIRECTA
directadental.com

If FenderMate® D-coat is separated during use, remove the matrix using a tweezer to avoid scratching gingiva.

602990	FenderMate D-Coat, refill, dark green, 18 pcs
602991	FenderMate D-Coat, refill, light green, 18 pcs
602992	FenderMate D-Coat, refill, dark blue, 18 pcs
602993	FenderMate D-Coat, refill, light blue, 18 pcs
602994	FenderMate D-Coat, assortment, 4 x 18 pcs

REF

Rx Only



106770-2019-03

FenderMate® (Registered Design and Patent pending) is a registered trademark owned by Directa AB.



DIRECTA AB

P.O. Box 723, Finvids väg 8
SE-194 27 Upplands Väsby,
Sweden

Tel: +46 8 506 505 75

DIRECTA INC

64 Barnabas Road, Unit 3
Newtown, CT 06470,
United States

203.491.2273