





FenderMate®

Matrix for Class II composite fillings

DIRECTA
directadental.com



FenderMate®

left regular

left narrow

right regular

right narrow



t= 0,05mm

Food approved plastic (Polypropene) and colours. Stainless steel.

If FenderMate® is separated during use, remove the matrix using a tweezer to avoid scratching gingiva.

602804	FenderMate Value Pack, dark green, 100 pcs
602805	FenderMate Value Pack, light green, 100 pcs
602806	FenderMate Value Pack, dark blue, 100 pcs
602807	FenderMate Value Pack, light blue, 100 pcs
602771	FenderMate, refill, dark green, 18 pcs
602772	FenderMate, refill, light green, 18 pcs
602773	FenderMate, refill, dark blue, 18 pcs
602774	FenderMate, refill, light blue, 18 pcs
602770	FenderMate, assortment, 4 x 18 pcs

Rx Only



FenderMate® (Registered Design and Patent pending) is a registered trademark owned by Directa AB.



DIRECTA AB

P.O. Box 723, Finvids väg 8
SE-194 27 Upplands Väsby,
Sweden

Tel: +46 8 506 505 75

DIRECTA INC

64 Barnabas Road, Unit 3
Newtown, CT 06470,
United States

203.491.2273

106854-2019-04-12