

DIRECTA CONTACTEZ



EN

ContacEZ Strip Restorative System

With three-steps - Adjust - Clean - Polish - ContacEZ Strip Restorative System provides an easy and effective way to complete crown seating. It breaks the process down into three simple steps: Adjustment, cleaning and polishing. Autoclavable.

1 Adjust proximal contact

The Black Diamond Dental Strip – The Proximal Contact Adjuster – Fine Diamond – Single-sided

- Adjust proximal contacts between the indirect restoration and the adjacent teeth.
- A strip to be used prior to definitive cementation of crowns, inlays/onlays, and veneers.
- Adjust and contour proximal contacts on composite fillings.

The Orange Serrated Diamond Strip

- A strip to be used to adjust and contour proximal contacts for composite fillings.
- The serrated edge readily separates and penetrates tight contacts to prepare for easy access to proximal contact adjustments.



REF	Product	Quantity
31108	Restorative System - Black	8 pcs
31132	Restorative System - Black	32 pcs



REF	Product	Quantity
31208	Restorative System - Orange	8 pcs
31232	Restorative System - Orange	32 pcs



REF	Product	Quantity
31408	Restorative System - White	8 pcs
31432	Restorative System - White	32 pcs



REF	Product	Quantity
31508	Restorative System - Blue	8 pcs
31532	Restorative System - Blue	32 pcs



REF	Product	Quantity
31308	Restorative System - Grey	8 pcs
31332	Restorative System - Grey	32 pcs

REF	Product	Content	Quantity
31008	Restorative System Assorted	2 Black, 1 Orange, 1 Gray, 2 White, 2 Blue	8 pcs
31040	Restorative System Assorted	8 Black, 8 Orange, 8 Gray, 8 White, 8 Blue	40 pcs
31023	Signature Set	3 Black, 3 Orange, 2 Gray, 3 White, 3 Blue, 2 Gold, 2 Copper, 2 Blue Narrow, 2 Pink	22 pcs

2 Clean out trapped cement

The White Serrated Strip – The Gentle Saw - Interproximal Cement Cleaner – No Diamond – Double-serrated.

- An ultra-thin serrated strip that is designed to safely cut and clean out excess cement.
- A strip to be used for interproximal cleaning of multiple veneers or composite fillings.

The Blue Serrated Strip – The Heavy-Duty Dental Saw – No Diamond – Double-serrated

- This heavy-duty dental saw strip is designed for more rigorous cutting and cleaning in the interproximal space.
- A strip that safely cuts into massive excess resin cement trapped in the interproximal space.

3 Polish to a natural finish

The Grey Final Polishing Strip – The Proximal Surface Polisher – Extra-fine Diamond – Single-sided

- A strip to be used to polish proximal contact surfaces of restorations and to restore a final finish in one step.

ContacEZ Strip IPR System


The **ContacEZ Strip IPR System** is an innovative set of dental instruments specifically crafted to make interproximal reduction safe, fast and easy. These flexible dental strips are ergonomically designed, patient-friendly, and conveniently labelled with specifications on every strip. ContacEZ’s interproximal reduction strips make IPR safe and accurate without creating sharp corners or cutting lips and gums. Designed for use with clear aligners and orthodontic treatments, these strips prevent excessive enamel removal and patient comfort and safety. ContacEZ IPR strips are autoclavable.

ContacEZ Rainbow Set

The **ContacEZ Rainbow Set** includes both the IPR basic strips as well as IPR PLUS strips. This 24-strip set gives dentists a full range of options for every possible case they encounter. Designed for use with clear aligners and other orthodontic treatments, these diamond strips are designed to be flexible, preserving the natural contours of the teeth while also preventing any harm or discomfort to patients.

REF	Product	Content	Quantity
32024	IPR Rainbow set	3 Yellow, 3 Red, 3 Blue, 3 Green, 2 Clear, 2 Cyan, 2 Purple, 2 Brown, 1 Mint, 1 Aqua, 1 Turquoise, 1 Teal, (0.07 - 0.50mm)	24 pcs

- **Solutions for most cases** – With the ContacEZ Strip IPR System, dentist have precise, easy-to-use devices to fit most patients’ IPR need.
- **Well contoured contacts** – The good-quality, stainless steel strips are flexible to preserve and conform to the natural contours of teeth during IPR. They are also taut enough to effortlessly access even tight contacts.
- **Accurate** – Paired with the Incremental Thickness Gauges, ContacEZ provides systematic, precise enamel removal without cutting lips, gums or cheeks.



IPR Strip System



Yellow Starter

0.07 mm Single-Sided.
Extra-Fine Diamond. Serrated

REF	Product	Quantity
32108	IPR Yellow	8 pcs
32132	IPR Yellow	32 pcs



Blue Widener

0.15 mm Double-Sided.
Coarse Diamond


REF	Product	Quantity
32508	IPR Blue	8 pcs
32532	IPR Blue	32 pcs



Red Opener

0.12 mm Double-Sided.
Medium Diamond

REF	Product	Quantity
32408	IPR Red	8 pcs
32432	IPR Red	32 pcs



Green Extra Widener

0.20 mm Double-Sided.
Coarse Diamond

REF	Product	Quantity
32608	IPR Green	8 pcs
32632	IPR Green	32 pcs

IPR PLUS



Clear Opener

0.10 mm Single-Sided.
Medium-Fine Diamond. Serrated

REF	Product	Quantity
32208	IPR PLUS Clear	8 pcs
32232	IPR PLUS Clear	32 pcs



Cyan Opener

0.15 mm Single-Sided.
Medium-Coarse Diamond

REF	Product	Quantity
32308	IPR PLUS Cyan	8 pcs
32332	IPR PLUS Cyan	32 pcs



Purple Widener

0.25 mm Double-Sided.
Coarse Diamond

REF	Product	Quantity
32708	IPR PLUS Purple	8 pcs
32732	IPR PLUS Purple	32 pcs



Brown Widener

0.30 mm Double-Sided.
Coarse Diamond

REF	Product	Quantity
32808	IPR PLUS Brown	8 pcs
32832	IPR PLUS Brown	32 pcs



Mint Widener

0.35 mm Double-Sided.
Coarse Diamond

REF	Product	Quantity
36608	IPR PLUS Mint	8 pcs
36632	IPR PLUS Mint	32 pcs



Aqua Widener

0.40 mm Double-Sided.
Coarse Diamond

REF	Product	Quantity
36708	IPR PLUS Aqua	8 pcs
36732	IPR PLUS Aqua	32 pcs



Turquoise Widener

0.45 mm Double-Sided.
Coarse Diamond

REF	Product	Quantity
36808	IPR PLUS Turquoise	8 pcs
36832	IPR PLUS Turquoise	32 pcs



Teal Widener

0.50 mm Double-Sided.
Coarse Diamond

REF	Product	Quantity
36908	IPR PLUS Teal	8 pcs
36932	IPR PLUS Teal	32 pcs

IPR PLUS Assorted Kits

REF	Product	Content	Quantity
32081	IPR PLUS Assorted Kit (0.10 mm - 0.50 mm) Small	1 Clear, 1 Cyan, 1 Purple, 1 Brown, 1 Mint, 1 Aqua, 1 Turquoise, 1 Teal	8 pcs
32082	IPR PLUS Assorted Kit (0.10 mm - 0.50 mm) Medium	2 Clear, 2 Cyan, 2 Purple, 2 Brown, 2 Mint, 2 Aqua, 2 Turquoise, 2 Teal	16 pcs
32084	IPR PLUS Assorted Kit (0.10 mm - 0.50 mm) Large	4 Clear, 4 Cyan, 4 Purple, 4 Brown, 4 Mint, 4 Aqua, 4 Turquoise, 4 Teal	32 pcs
32009	IPR PLUS Optional Assorted Kit (0.10 mm - 0.25 mm)	2 Clear, 2 Cyan, 2 Purple, 2 Brown	8 pcs
36008	IPR PLUS Optional Assorted Kit (0.35 mm - 0.50 mm)	2 Mint, 2 Aqua, 2 Turquoise, 2 Teal	8 pcs

ContacEZ Strip Subgingival System

The Subgingival strips are designed to remove carefully and conveniently hard-to-reach subgingival rough margins and overhangs and polish the surface to restore a natural finish. This system is specially designed to access the subgingival area to make overhang removal simple while preventing damage to patients’ gingival tissue and lips. In addition, these ultra-narrow strips are perfect for cleaning out cement debris in the subgingival area of the interproximal space after crown cementation. ContacEZ Subgingival strips are autoclavable.



REF	Product	Content	Quantity
33108	Subgingival Strip Assorted	2 Blue, 3 Gold, 3 Copper	8 pcs
33124	Subgingival Strip Assorted	8 Blue, 8 Gold, 8 Copper	24 pcs



REF	Product	Quantity
33708	Subgingival Strip BLUE	8 pcs
33732	Subgingival Strip BLUE	32 pcs



REF	Product	Quantity
33408	Subgingival Strip GOLD	8 pcs
33432	Subgingival Strip GOLD	32 pcs



REF	Product	Quantity
33608	Subgingival Strip COPPER	8 pcs
33632	Subgingival Strip COPPER	32 pcs

Cut

The Blue Narrow Strip – 0.07 mm – Subgingival Cutter Double-serrated – No Diamond

- The Blue Narrow Strip can be used when large, protruding overhangs are detected.
- The thin, serrated strip easily accesses the subgingival area and cuts away excess composite material or crown cement.

Trim

The Gold Narrow Strip – 0.09 mm – Subgingival Trimmer - Single- sided – Medium Coarse Diamond – Central Gateway

- A strip to be used when overhangs are smaller or do not stick out from the tooth.
- The thicker diamond strip allows you to quickly trim away overhangs, while the blank central gateway enables smooth access to the interproximal space.

Polish

The Copper Narrow Strip – 0.06 mm – Subgingival Polisher Single-sided – Extra Fine Diamond

- A strip to be used to clean and polish the restoration after removing overhang.
- Restores a natural finish to the surface of the tooth in one easy step.

ContacEZ Strip Hygienist System

Remove stains, discolorations, and clean proximal contact areas without damaging enamel surfaces.

The ContacEZ Hygienist strip is a perforated 0.05 mm thin dental strip designed to remove extrinsic discoloration and plaque from interproximal surfaces without causing damage to tooth enamel. The non-abrasive, flexible, stainless-steel strips clean proximal contact areas of tough stains from coffee, tea, tobacco, and red wine. ContacEZ Hygienist strips are autoclavable.



REF	Product	Quantity
35108	Hygienist Strip	8 pcs
35132	Hygienist Strip	32 pcs

Strip Tray

These convenient, autoclavable Strip trays are perfect for keeping your strips organized and accessible. Each tray holds 8 strips.

REF	Product	Quantity
92501	Strip Tray	1 pc
92504	Strip Tray	4 pcs



Incremental thickness gauge



The **ContacEZ Incremental Thickness Gauge** efficiently and accurately confirms the amount of interproximal enamel removal or contact reduction for orthodontic treatment.

REF	Product	Quantity
92106	Incremental Thickness Gauge	1 pc

Design By Dentists

Directa's products are developed, tested and evaluated by our selected team of highly-qualified and renowned Swedish and international dental professionals. The aim is to make life easier for all dental professionals in their daily clinical work.

DIRECTA



DIRECTA AB

Finvids väg 8
194 47 Upplands Väsby, SWEDEN
+46 8 506 505 75
info@directadental.com
www.directadental.com

DIRECTA INC

64 Barnabas Road, Unit 3
Newtown, CT. 06470, USA
+1 203.491.2273
infousa@directadental.com
www.directausa.com



directadental



directadental



DirectaAB

The content of this brochure, including layout and images, is protected by copyright. ©2026 Directa AB. All rights reserved.

