

DIRECTA
PractiPal®





What are the ideal features of a tray system?

Infection control is crucial, especially when it comes to instrument reprocessing. To ensure efficiency in cleaning, disinfection, and sterilization, you need a good tray system that can handle instruments effectively. A good tray system should allow for efficient handling of instruments throughout the reprocessing/cleaning cycle.

A tray system with pre-prepared instrument sets used for each different type of treatment prevents the risk of confusing clean and used instruments, saves time, and improves the clinic's efficiency. Instruments should be held securely within the tray, but also free to move for thorough cleaning. Trays that hold instruments too firmly or

have insufficient water penetration can leave microbiological residue on the instruments.

Logistics of instrument handling are often overlooked but are important. Overloading can make it difficult to keep track of instruments. PractiPal solves this issue with its clever tray design, making it easier for you to stay organized and maintain high standards of infection control. Regular reviews are recommended to ensure cleaning efficiency, storage capacity, and sterility are optimized. Overall, PractiPal is an ideal tray system that can optimize your infection control process.



**Experience the
convenience and
effectiveness with
PractiPal® in your
dental clinic.**

This lightweight and colourful tray is both a pleasure to use and exceptionally durable, thanks to its material, which does not absorb heat in autoclaves or washer disinfectors. This also saves energy costs and aids in the cleaning process. PractiPal Tray System includes trays with bur stands, consumables, and accessories that can easily be secured onto the trays, making the entire cleaning and disinfection process more efficient.

The tray clamp and bur stand are both easy to open, and the fibreglass-reinforced plastic ensures that instruments are held securely, reducing

the risk of accidents, and simplifying instrument inspection.

PractiPal allows for the efficient disposal of waste and holds instruments and burs in place for cleaning and sterilization. By using pre-prepared trays containing only what is necessary for the procedure, you can easily cycle through tray configurations. This eliminates the need to search for bur and accessory holders that can easily be misplaced. Once treatment is complete, the used tray can be replaced with a new configuration for the next patient.

Organized work flow

Pre-configured tray systems minimize the quantity of instruments required whilst maximizing treatment efficiency, thanks to a simple cleaning procedure.

Safety first

PractiPal holds all instruments securely, preventing accidents and cross-contamination.

Improved hygiene control

The open design allows for maximum water and steam penetration, ensuring an optimal cleaning effect. Instruments are secure yet loose enough to avoid contact areas between the instruments and tray.

Light and non-metallic

Avoids electrolytic corrosion. Lightness of the tray reduces heat energy loss and risk of overload in the autoclave.



PractiPal® do not compromise on predictable cleaning.

The belief that sterile is “clean and safe” is a common misconception. Residue can still remain on instruments even after cleaning, and this “sterile” residue can cause severe inflammatory reactions. Predictable cleaning is essential for 100% control. PractiPal trays, with a minimized frame, incorporate an ingenious means of securing instruments loosely. This ensures maximum water penetration for optimal washing, ultrasonic cleaning or simply rinsing.

Optimal cleaning, easy storage and safe handling makes PractiPal the product of choice. PractiPal includes bur, file and accessory holders that clip securely into the tray but are also easy to remove. The wide range of colours and accessories make it easy for you to design unique configurations to suit each surgery, dentist and individual operator. Once a cassette is configured, it can be reused time and time again.



PractiPal® Full Tray

PractiPal® Full Tray

(280x188x37 mm)

Two half pieces assembled with joining clips. Organises and secures storage of instruments. Holds 9 instruments. Bur/Endo stand holds 16 burs or 28 files. Available in nine colours. Colour of the tray: black. PractiPal Complete Set contains: Tray (2 pcs), clips (4 pcs), Instrument Clamp (1 pc), Bur/Endo stand for 16 burs or 28 files (1 pc), MultiBlock (1 pc), MultiDappen (5 pcs), MultiPoint (5 pcs), Waste Cup (5 pcs), Foam Large (1 pc) and Foam Small (5 pcs).

Prod No	Product Description
115199	PractiPal Complete Set, black/grey
115191	PractiPal Complete Set, black/red
115192	PractiPal Complete Set, black/yellow
115193	PractiPal Complete Set, black/blue
115194	PractiPal Complete Set, black/green
115195	PractiPal Complete Set, black/lilac
115196	PractiPal Complete Set, black/pink
115197	PractiPal Complete Set, black/mint green
115198	PractiPal Complete Set, black/orange



PractiPal® Half Tray

PractiPal® Half Tray

(140x188x37 mm)

PractiPal Half Tray suits all set-ups and contains one half tray and one instrument clamp, holding up to 9 instruments. Available in nine colours. Colour of the tray: black

Prod No	Product Description
115100	PractiPal Tray, 2 clips
115106	PractiPal Tray and Instrument Clamp, black/lilac
115107	PractiPal Tray and Instrument Clamp, black/pink
115108	PractiPal Tray and Instrument Clamp, black/mint green
115109	PractiPal Tray and Instrument Clamp, black/orange
115110	PractiPal Tray and Instrument Clamp, black/grey
115111	PractiPal Tray and Instrument Clamp, black/red
115112	PractiPal Tray and Instrument Clamp, black/yellow
115113	PractiPal Tray and Instrument Clamp, black/blue
115114	PractiPal Tray and Instrument Clamp, black/green



PractiPal® Mini Tray


PractiPal® Mini Tray

(80x188x37 mm)

PractiPal Mini Tray holds 5 instruments. This makes it ideal for check-ups, for dental hygienists and to stock spare instruments. Available in nine colours. Colour of the tray: black.

Prod No	Product Description
115020	PractiPal Mini Tray, Black
115050	PractiPal Mini Tray and Mini Clamp, black/grey
115051	PractiPal Mini Tray and Mini Clamp, black/red
115052	PractiPal Mini Tray and Mini Clamp, black/yellow
115053	PractiPal Mini Tray and Mini Clamp, black/blue
115054	PractiPal Mini Tray and Mini Clamp, black/green
115055	PractiPal Mini Tray and Mini Clamp, black/lilac
115056	PractiPal Mini Tray and Mini Clamp, black/pink
115057	PractiPal Mini Tray and Mini Clamp, black/mint green
115058	PractiPal Mini Tray and Mini Clamp, black/orange

Select your
colour coded
treatment tray



PractiPal® Instrument Clamp

Secures instruments for optimal cleaning without damage. Available in nine different colours and two sizes. One for Full and Half Tray and one for PractiPal Mini Tray.

Prod No	Product Description
115150	PractiPal Instrument Clamp, grey
115151	PractiPal Instrument Clamp, red
115152	PractiPal Instrument Clamp, yellow
115153	PractiPal Instrument Clamp, blue
115154	PractiPal Instrument Clamp, green
115155	PractiPal Instrument Clamp, lilac
115156	PractiPal Instrument Clamp, pink
115157	PractiPal Instrument Clamp, mint green
115158	PractiPal Instrument Clamp, orange

Prod No	Product Description
115040	PractiPal Mini Clamp, grey
115041	PractiPal Mini Clamp, red
115042	PractiPal Mini Clamp, yellow
115043	PractiPal Mini Clamp, blue
115044	PractiPal Mini Clamp, green
115045	PractiPal Mini Clamp, lilac
115046	PractiPal Mini Clamp, pink
115047	PractiPal Mini Clamp, mint green
115048	PractiPal Mini Clamp, orange



PractiPal® Compact Bur Stand

Practical and safe bur holder for organization of burs.
Holds 2 long burs, 5 RA and 8 FG burs.
Easy handling - Minimizes sharps injuries.
Can be used together with the PractiPal system or as a stand alone.

Prod No: 115121 PractiPal Compact Bur Stand, lilac, 1 pc
 Prod No: 115122 PractiPal Compact Bur Stand, pink, 1 pc
 Prod No: 115123 PractiPal Compact Bur Stand, mint green, 1 pc
 Prod No: 115124 PractiPal Compact Bur Stand, orange, 1 pc
 Prod No: 115185 PractiPal Compact Bur Stand, grey, 1 pc
 Prod No: 115186 PractiPal Compact Bur Stand, red, 1 pc
 Prod No: 115187 PractiPal Compact Bur Stand, yellow, 1 pc
 Prod No: 115188 PractiPal Compact Bur Stand blue, 1 pc
 Prod No: 115189 PractiPal Compact Bur Stand, green, 1 pc



PractiPal® Compact File Stand

Compact file holder for the practical and safe organization of endodontic files.
The file stand holds 13 files.
Easy handling - Minimizes sharps injuries.
Can be used together with the PractiPal system or as a stand alone.

Prod No: 115176 PractiPal Compact File Stand, lilac, 1 pc
 Prod No: 115177 PractiPal Compact File Stand, pink, 1 pc
 Prod No: 115178 PractiPal Compact File Stand, mint green, 1 pc
 Prod No: 115179 PractiPal Compact File Stand, orange, 1 pc
 Prod No: 115180 PractiPal Compact File Stand, grey, 1 pc
 Prod No: 115181 PractiPal Compact File Stand, red, 1 pc
 Prod No: 115182 PractiPal Compact File Stand, yellow, 1 pc
 Prod No: 115183 PractiPal Compact File Stand, blue, 1 pc
 Prod No: 115184 PractiPal Compact File Stand, green, 1 pc



PractiPal® Bur/Endostand

Functioning as a bur stand or file holder, the same stand can hold either 28 files or 16 burs. Allows for easy inspection of files and burs, minimizing the risk of sharps injury.

Prod No	Product Description
115130	PractiPal Bur/Endostand, grey
115131	PractiPal Bur/Endostand, red
115132	PractiPal Bur/Endostand, yellow
115133	PractiPal Bur/Endostand, blue
115134	PractiPal Bur/Endostand, green
115135	PractiPal Bur/Endostand, lilac
115136	PractiPal Bur/Endostand, pink
115137	PractiPal Bur/Endostand, mint green
115138	PractiPal Bur/Endostand, orange



PractiPal® MultiBlock

Holder that stores disposable items within its multiple inserts. Fits PractiPal Half and Full Tray.

PractiPal MultiBlock contain: PractiPal MultiBlock (1 pc), PractiPal MultiDappen (5 pcs), PractiPal MultiPoint (5 pcs), PractiPal Waste Cup (5 pcs) and PractiPal Foam small (5 pcs).

Prod No	Product Description
115140	PractiPal MultiBlock, grey
115141	PractiPal MultiBlock, red
115142	PractiPal MultiBlock, yellow
115143	PractiPal MultiBlock, blue
115144	PractiPal MultiBlock, green
115145	PractiPal MultiBlock, lilac
115146	PractiPal MultiBlock, pink
115147	PractiPal MultiBlock, mint green
115148	PractiPal MultiBlock, orange

PractiPal® Half Tray Compact File Set

(140x188x37 mm)

Compact file holder stores 13 files. Capacity: 5 instruments. Available in nine colours. Colour of the tray: black. PractiPal Compact File Set contains: PractiPal Half Tray, PractiPal Instrument Clamp for Half Tray, Compact File Stand and Waste Cups (4 pcs).

Prod No	Product Description
115200	PractiPal Compact File Set, black/grey
115201	PractiPal Compact File Set, black/red
115202	PractiPal Compact File Set, black/yellow
115203	PractiPal Compact File Set, black/blue
115204	PractiPal Compact File Set, black/green
115205	PractiPal Compact File Set, black/lilac
115206	PractiPal Compact File Set, black/pink
115207	PractiPal Compact File Set, black/mint green
115208	PractiPal Compact File Set, black/orange



PractiPal® Half Tray Compact Bur Set

(140x188x37 mm)

Compact Bur Set holds 2 extra long burs, 5HP/RA and 8FG. Capacity: 5 instruments. Available in nine colours. Colour of the tray: black. PractiPal Compact Bur Set contains: PractiPal Half Tray, PractiPal Instrument Clamp for Half Tray, Compact Bur Stand and Waste Cups (4 pcs).

Prod No	Product Description
115160	PractiPal Compact Bur Set, black/red
115161	PractiPal Compact Bur Set, black/yellow
115162	PractiPal Compact Bur Set, black/green
115163	PractiPal Compact Bur Set, black/blue
115164	PractiPal Compact Bur Set, black/grey
115170	PractiPal Compact Bur Set, black/lilac
115171	PractiPal Compact Bur Set, black/pink
115172	PractiPal Compact Bur Set, black/mint green
115173	PractiPal Compact Bur Set, black/orange





PractiPal MultiDappen, MultiPoint and Waste Cup contain a proportion of renewable wood fibres. The natural origin of the fibres give them a light wood-like scent, which is completely normal.

PractiPal® Waste Cup

Practical container for debris or liquids. Instruments may be scraped against the special jagged edges to remove debris. Can be used stand-alone or placed in the MultiBlock.

Prod No	Product Description
115002	PractiPal Waste Cup, 300 pcs

PractiPal® MultiPoint

Storage of Paper Points, Guttapercha Points and fluid. Can be used alone or placed in the MultiBlock.

Prod No	Product Description
115001	PractiPal MultiPoint, 300 pcs

PractiPal® MultiDappen

Designed with indentations for different volumes of bonding agents, primers, etc. Can be used alone or placed on the MultiBlock.

Prod No	Product Description
115000	PractiPal MultiDappen, 300 pcs



PractiPal® Clips

Used to assemble two PractiPal Half Tray pieces together.

Prod No	Product Description
115120	PractiPal Clips, 2 pcs



PractiPal® Foam

The unique slot design of the foam scraper allows for easy removal of debris from files and instruments. Large: 50x32x23 mm. Small: 30x15x13 mm.

Prod No	Product Description
115003	PractiPal Foam, small, 30x15x13 mm, 300 pcs
115004	PractiPal Foam, large, 50x32x23 mm, 50 pcs

Design By Dentists

Directa's products are developed, tested and evaluated by our selected team of highly-qualified and renowned Swedish and international dental professionals. The aim is to make life easier for all dental professionals in their daily clinical work.

DIRECTA



DIRECTA AB

Finvids väg 8
194 47 Upplands Väsby, SWEDEN
+46 8 506 505 75
info@directadental.com
www.directadental.com

DIRECTA INC

64 Barnabas Road, Unit 3
Newtown, CT. 06470, USA
+1 203.491.2273
infousa@directadental.com
www.directausa.com

 directadental

 directadental

 DirectaAB

The content of this brochure, including layout and images, is protected by copyright. ©2026 Directa AB. All rights reserved.

