

ContacEZ Strips IPR System

Solutions for most cases – The **ContacEZ Strips IPR System** provides the dental professional with accurate, easy-to-use devices to meet most patient IPR needs.

Contoured contacts - The high quality stainless steel strips are flexible to preserve and conform to the natural contours of the teeth during IPR. They are also tight enough to reach even the tightest contacts with ease.

Accurate – Paired with the Incremental Thickness Gauges, ContacEZ provides systematic, precise enamel removal without cutting lips, gums or cheeks.



The **ContacEZ Strips IPR System** is an innovative set of dental instruments specifically crafted to make interproximal reduction safe, fast, and easy. These flexible dental strips are ergonomically designed, patient-friendly and each strip is marked with the thickness of the strip and whether it is single-sided (S) or double-sided (D). ContacEZ's interproximal reduction strips make IPR safer and more accurate without creating sharp corners or cutting lips and gums. Designed for use with clear aligners and orthodontic treatments, these strips prevent excessive enamel removal and maximise patient comfort and safety. Autoclavable and reusable.

ContacEZ Strip IPR System



Yellow Starter

0.07 mm
Single-Sided. Extra-Fine Diamond. Serrated
Prod No: 32108 8 pcs, 32132 32 pcs



Red Opener

0.12 mm
Double-Sided. Medium Diamond
Prod No: 32408 8 pcs, 32432 32 pcs



Blue Widener

0.15 mm
Double-Sided. Coarse Diamond
Prod No: 32508 8 pcs, 32532 32 pcs



Green Extra-Widener

0.20 mm
Double-Sided. Coarse Diamond
Prod No: 32608 8 pcs, 32632 32 pcs

ContacEZ Strip IPR PLUS System



Clear Opener

0.10 mm
Single-Sided. Med-Fine Diamond. Serrated
Prod No: 32208 8 pcs, 32232 32 pcs



Cyan Opener

0.15 mm
Single-Sided. Medium Coarse Diamond
Prod No: 32308 8 pcs, 32332 32 pcs



Purple Widener

0.25 mm
Double-Sided. Coarse Diamond
Prod No: 32708 8 pcs, 32732 32 pcs



Brown Widener

0.30 mm
Double-Sided. Coarse Diamond
Prod No: 32808 8 pcs, 32832 32 pcs



Mint Widener

0.35 mm
Double-Sided. Coarse Diamond
Prod No: 36608 8 pcs, 36632 32 pcs



Aqua Widener

0.40 mm
Double-Sided. Coarse Diamond
Prod No: 32708 8 pcs, 32732 32 pcs



Turquoise Widener

0.45 mm
Double-Sided. Coarse Diamond
Prod No: 36808 8 pcs, 36832 32 pcs



Teal Widener

0.50 mm
Double-Sided. Coarse Diamond
Prod No: 32908 8 pcs, 32932 32 pcs


Prod No Product Description

32008 IPR Assorted (0.07mm - 0.20mm) Small, 2 Yellow, 2 Red, 2 Dark Blue, 2 Green. 8 pcs

32024 IPR Rainbow Set (0.07mm - 0.50mm) 3 Yellow, 3 Red, 3 Blue, 3 Green, 2 Clear, 2 Cyan, 2 Purple, 2 Brown, 1 Mint, 1 Aqua, 1 Turquoise, 1 Teal. 24 pcs

ContacEZ Strips Restorative System

The three-step **ContacEZ Strips Restorative System** offers you an easy and effective way to finish a crown seating, breaking the process down into three simple steps: **Adjust**, **Clean**, and **Polish**. Autoclavable and reusable.



Prod No	Product Description
31008	ContacEZ Strips Restorative System Assorted. 8 pcs 2 Black, 1 Orange, 1 Grey, 2 White, 2 Blue
31040	ContacEZ Strips Restorative System Assorted. 40 pcs 8 Black, 8 Orange, 8 Grey, 8 White, 8 Blue
31023	ContacEZ Strips Restorative Signature Set. 22 pcs 3 Black, 3 Orange, 2 Grey, 3 White, 3 Blue, 2 Gold, 2 Copper, 2 Blue Narrow, 2 Pink

Adjust



Black - Diamond Strip
Single-Sided. Fine Diamond
Prod No: 31108 8 pcs, 31132 32 pcs



Orange - Serrated Diamond Strip
Single-Sided. Extra-Fine Diamond
Double Serrated
Prod No: 31208 8 pcs, 31232 32 pcs

Clean



White - Gentle Saw
Double Serrated. No Diamond
Prod No: 31408 8 pcs, 31432 32 pcs



Blue - Heavy Duty Saw
Double Serrated. No Diamond
Prod No: 31508 8 pcs, 31532 32 pcs

Polish



Grey - Final Polishing Strip
Single-Sided. Ultra-Fine Diamond
Prod No: 31308 8 pcs, 31332 32 pcs

ContacEZ Strips Subgingival System

The **ContacEZ Strips Subgingival System** is designed to carefully and conveniently remove hard-to-reach subgingival rough margins and overhangs. They provide easy access to the subgingival area while avoiding damaging patient's gingival tissue and lips. In addition, these ultra-narrow strips are ideal for cleaning out cement debris in the subgingival area of the interproximal space after crown cementation. Autoclavable and reusable.



Prod No	Product Description
33108	ContacEZ Strips Subgingival System Assorted. 8 pcs 2 Blue, 3 Gold, 3 Copper
33124	ContacEZ Strips Subgingival System Assorted. 24 pcs 8 Blue, 8 Gold, 8 Copper



Gold
Single-Sided. Medium Diamond
Prod No: 33408 8 pcs, 33432 32 pcs



Blue
No Diamond. Double Serrated
Prod No: 33708 8 pcs, 33732 32 pcs



Copper
Single-Sided. Extra-Fine Diamond
Prod No: 33608 8 pcs, 33632 32 pcs

ContacEZ Strips Hygienist System

Remove stains and clean proximal contact areas without damaging enamel surfaces. The **ContacEZ Hygienist Strip** is an advanced, high-precision dental strip designed to remove extrinsic discoloration and plaque from interproximal surfaces without causing damage to tooth enamel. The non-abrasive, flexible, stainless-steel strips clean proximal contact areas of tough stains from coffee, tea, tobacco, and red wine. It is quick, safe and ensure patient satisfaction and save hygienist chair time. Autoclavable and reusable.



Pink
Perforated. Stainless Steel
Prod No: 35108 8 pcs, 35132 32 pcs

35008 **ContacEZ Strips Hygienist Removal Kit**
2 Subgingival Blue, 2 Subgingival Copper,
2 Subgingival Gold, 2 Hygienist Strips. 8 pcs

Video: ContacEZ -
In the clinic with Cecilia Ömo

