

EASIER - BETTER - FASTER




TrollFoil takes the guesswork out of occlusal adjustments.

TrollFoil can be used under a wide variety of clinical situations including wet or dry teeth, limited opening, limited vestibular space, gaggers, and metal and non-metallic restorations. You are able to verify occlusal contacts for both arches to an entire quadrant. The double-sided foil is only 8 microns thick, and it has no problem marking wet surfaces, dry surfaces or highly polished surfaces such as cast gold or BruxZir.

No need for forceps
The foil is pre-mounted in
a sturdy plastic frame.

No ink on your fingers
The foil is covered on both
sides at the grip.




Protective strip
Easily removed before use.

Super thin foil
The foil is only 8 micron thick
and double-sided.
Available in blue or red.

When using TrollFoil there is no need for forceps as TrollFoil is pre-mounted in a sturdy plastic frame keeping the foil stretched and straight allowing you to hold the foil with your fingers without getting colour on them.

The double-sided foil is covered on both sides with a protective strip you easily remove before use.



Video: TrollFoil
"I have been using TrollFoil for 10 years..." - Folke Ullberg



Video: TrollFoil
The only articulating foil you will ever need!

Scan the QR-code with your QR-reader on your mobile phone to watch the videos.



TrollFoil Mini for
one or two teeth

Art.

12000402 TrollFoil Mini, 100 pcs



TrollFoil 4.5 - ultra thin
only 4.5 micron

Art.

12000145 TrollFoil 4.5, 80 pcs



TrollFoil - super thin
only 8 micron

Art.

12000501 TrollFoil - blue, 500 pcs
12000600 TrollFoil - red, 500 pcs
12000101 TrollFoil - blue, 100 pcs
12000700 TrollFoil - red, 100 pcs



For more information visit
www.directadental.com or
scan the QR-Code.

DIRECTA