


Directa Crowns Molars Mould Guide


Translucent/Translucide/Translucida/Translucenti

DIRECTA Directa Temporary Crown


Upper Right/Höger ÖK/Haut droit/OK-Rechts/Superior Derecha/Superiore Destro




Upper Left/Vänster ÖK/Haut gauche/OK-Links/Superior Izquierda/Superiore Sinistro



Lower Right/Höger UK/Bas droite/UK-Rechts/Inferior Derecha/Inferiore Destro



Lower Left/Vänster UK/Bas gauche/UK-Links/Inferior Izquierda/Inferio Sinistro



Opaque/Opak/Opaca/Opache
Translucent/Translucide/Translucida/Translucenti

Centrals/Centraler/Centrales/OK-Front/Centrales/Centrali

mm	14.0	13.0	12.4	12.0	11.5	11.0	10.6	10.6	11.0	11.5	12.0	12.4	13.0	14.0
mm	10.0	9.2	8.7	8.4	8.2	7.8	7.5	7.5	7.8	8.2	8.4	8.7	9.2	10.0

centrals

Laterals/Lateraler/Latérales/OK-Front/Laterales/Laterali

Teeth	2	20	21	22	23	24	25	26	27	28	29	200
Height (mm)	11.0	10.4	10.1	9.5	8.8	8.2	8.2	8.8	9.5	10.1	10.4	11.0
Width (mm)	7.3	6.8	6.5	6.2	5.9	5.5	5.5	5.9	6.2	6.5	6.8	7.3

Cuspids/Hörntänder/Canines/Eckzähne/Caninos/Cuspidi

mm ↑	13.0	12.4	11.6	11.4	11.1	10.8	10.5	10.5	10.8	11.1	11.4	11.6	12.4	13.0
mm ↔	9.0	8.5	8.1	7.9	7.7	7.5	7.3	7.3	7.5	7.7	7.9	8.1	8.5	9.0
	300	301	30	31	32	33	34	35	36	37	38	39	302	303
	<i>cuspid</i>													

**Upper bicuspids/ÖK-premolarer/Prémolaires supérieures/
OK Prämolaren/Premolares superiores/Bicuspidi superiori**

mm	9.8	9.4	8.9	8.6	8.4	7.5	8.1	8.5	9.0	9.5
mm	7.4	6.9	6.7	6.3	6.0	6.0	6.3	6.7	7.2	7.4
	40	41	42	43	44	45	46	47	48	49

bicusplids upper

**Lower bicuspid/UK-premolarer/Prémolaires inférieures/
UK Prämolaren/Premolares inferiores/Bicuspidi inferiori**

mm	9.8	9.5	9.2	8.8	8.2	7.1	7.3	8.0	8.1	8.4
mm	7.4	7.2	6.8	6.6	6.3	6.4	6.8	7.1	7.5	8.0

50 51 52 53 54 55 56 57 58 59

bicuspids lower

**Lower incisors/UK-incisiver/Incisives inférieures/
UK Front/Incisivos inferiores/Incisivi inferiori**

mm	12.0	11.5	11.0	10.5	10.0	8.0	8.5	9.1	9.6	10.3
mm	6.2	5.6	5.3	5.1	4.9	4.5	4.9	5.2	5.5	5.8

Incisors lower